

## Making Competent Cells

Need: 1-50 ml Falcon tube of  $\text{CaCl}_2$   
5-15 ml minicentrifuge tubes  
2-15 ml Falcon tube

### Night Before

Prepare Stock of DH5  $\alpha$  cells

1. Add 5 ml of LB broth into 15 ml Falcon tube
2. Add loop of cells from plate
3. Incubate in shaking incubator at  $37^\circ\text{C}$  O/N

### Next Day

Add 100  $\mu\text{l}$  DH5  $\alpha$  mixture + 9.9 ml LB broth in 15 ml Falcon tube  
↓ 3-4 hours in shaking incubator

[Add 500  $\mu\text{l}$  of vial into 1.7 ml minicentrifuge tubes] x 5

↓ centrifuge at 10000 RPM for 5 min  
-if no pellet at base raise RPM or add more culture

Remove Supernatant x 5

↓

[Add 500  $\mu\text{l}$   $\text{CaCl}_2$  - resuspend solution] x 5 → 15 min on ice

↓ centrifuge at 10000 RPM for 5 min

Remove supernatant x 5

↓

[Add 100  $\mu\text{l}$   $\text{CaCl}_2$  - resuspend solution] x 5

↓

30 min on ice → add DNA within this time

Result 5 minicentrifuge tubes containing 100  $\mu\text{l}$  competent DH5  $\alpha$  cells

Transformation

- ① Made 100ng/ml original stock
- ② Diluted 1:10

Need: 5-1.7 mL vials + 3 vials chitinase enzymes made according to instructions + RFP (J04460) + J04500  
ON ICE ON ICE ON ICE



Heat Shock in floater at 42°C for 60 seconds

↓  
Ice in floater 5min

↓  
Add 500ul per tube LB broth

↓  
Incubate at 37°C for 60 min

↓  
Spin down for 30 seconds at 10000 RPM

↓  
Remove supernatant

↓  
Resuspend in 200ul LB broth

↓  
Plate 200ul of each vial in appropriate antibiotic plate → Incubate at 37°C O/N

J04500/J04450-p681c3  
-chlor resistant

Chi enzymes - pUC57  
-amp resistant

# Mini Prep

1-2 days  
before  
Mini Prep

Select single colony from freshly streaked transformed plate using inoculating loop

add to 10-15 mL LB media  
+  
5  $\mu$ L Amp in 50 mL Falcon tube

Incubate 12-16 hours at 37°C  
in Shaking incubator

## Prepare Reagents

1 mL of culture to 1.7 mL  
minicentrifuge tube

Spin at 8000 rpm  
for 1 min

Decant

Add 1 more mL of culture  
to 1.7 mL centrifuge tube

Spin at 8000 rpm  
for 10 min

Remove supernatant  
in beaker with bleach

Add 600  $\mu$ L solution  
III (neutralization buffer)  
invert until white precipitate  
forms

invert + rotate tube  
several times (get clear  
lysate)

Add 400  $\mu$ L  
solution II (lysis buffer)

thoroughly resuspend  
cells in 500  $\mu$ L  
RNase A

Spin at 14000 rpm  
for 10 min. Want to see  
pellet debris.




Insert HiBind DNA  
mini Column into 2 mL  
Collection tube (label all tubes)

Add 700  $\mu$ L clear lysate  
from spun tube aspirating  
(swirling around edge of mini  
column) into HiBind DNA  
mini column. Leave pellet undisturbed.

Spin at 14000 rpm for 1 min

Pull out mini column tubes, discard  
liquid, put column tubes back into  
collection tubes.

Next Page



### Dry spin

HiBind DNA minitubes at 14000 rpm for 3 min to dry the column. (no ethanol left)

\* Repeat 3 times to remove All wash buffer

Transfer HiBind DNA mini column to clean 1.5 mL microcentrifuge tube

Add 80-100  $\mu$ L Elution buffer to Centre of column membrane

Allow column to sit at room temperature for 1 min.

Spin at 14000 rpm for 2 min


↓  
This represents 70% bound DNA

# Gel Electrophoresis

Make gel agarose

- 1g agarose in 100 mL TAE
- Microwave 30s x 2 + 3<sup>rd</sup> time for 10s  
↳ saran wrap top with a hole in the top
- Wait until cool enough to touch for 3s
- Add 5  $\mu$ l SYBR green
- Swirl until colour is fully dissolved
- Pour into mould (covered by masking tape)
- Insert comb at end slot

- After gel is solidified (~30 mins), remove comb & tape
- Place into gel electrophoresis machine
- Fill with TAE until top of gel is covered



- Pipette all into each well (spin down first if necessary)
- 125 V for 40 min, checking every 10 min
- Ladder  $\rightarrow$  compared to  
~ 3700 bp for chitinases + PUC 57  
~ 3100 bp for J04500 + pSB1C3  
~ 2200 bp for J04450 + pSB1C3

# Restriction Digest

In 1.7 mL centrifuge tube, add:

30.5  $\mu$ l ddH<sub>2</sub>O  
5  $\mu$ l buffer cutsmart  
12.5  $\mu$ l chitinase inserts  
1  $\mu$ l Spe I  
1  $\mu$ l Pst I  
    ↙    ↓    ↘  
Pc1-1 P<sub>g</sub>1-1 P<sub>g</sub>1-2

In 1.7 mL centrifuge tube, add:

30.5  $\mu$ l ddH<sub>2</sub>O  
50  $\mu$ l buffer cutsmart  
12.5  $\mu$ l J04... vectors  
    J04450    {    J04500  
1  $\mu$ l XbaI    {    1  $\mu$ l Spe I  
1  $\mu$ l Pst I    {    1  $\mu$ l Pst I  
    ↓                    ↓  
    J04460                J04500

Pipette gently 5x

Centrifuge @ max for 10s

Incubate @ 37°C for 1 hr

remainder freeze @ -20°C for ligation


Gel electrophoresis to check if digest worked (want 2 bands - one of enzyme + one of backbone)

# Ligation

(F)

80°C for 20 min  
water bath to  
deactivate restriction  
enzymes.

J04450 J04500 PgeI-1 Pcl-1 PgeI-2




in each 0.6 mL tube add:

2  $\mu$ L vector (J04....)  
6.5  $\mu$ L insert (chitinases)  
1  $\mu$ L T4 DNA buffer  
0.5  $\mu$ L T4 DNA ligase  

---

10  $\mu$ L



leave to grow for 24 hours.