

Mini Prep

1-2 days
before
Mini Prep

Select single colony from freshly streaked transformed plate using inoculating loop

add to 10-15 mL LB media
+
5 μ L Amp in 50 mL Falcon tube

Incubate 12-16 hours at 37°C
in Shaking incubator

Prepare Reagents

1 mL of culture to 1.7 mL
minicentrifuge tube

Spin at 8000 rpm
for 1 min

Decant

Add 1 more mL of culture
to 1.7 mL centrifuge tube

Spin at 8000 rpm
for 10 min

Remove supernatant
in beaker with bleach

Add 600 μ L solution
III (neutralization buffer)
invert until white precipitate
forms

invert + rotate tube
several times (get clear
lysate)

Add 400 μ L
solution II (lysis buffer)

thoroughly resuspend
cells in 500 μ L
RNase A

Spin at 14000 rpm
for 10 min. Want to see
pellet debris.

Insert HiBind DNA
mini Column into 2 mL
Collection tube (label all tubes)

Add 700 μ L clear lysate
from spun tube aspirating
(swirling around edge of mini
column) into HiBind DNA
mini column. Leave pellet undisturbed.

Spin at 14000 rpm for 1 min

Pull out mini column tubes, discard
liquid, put column tubes back into
collection tubes.

Next Page