

Extract Genomic DNA (gDNA) from '288c' 2/27/14

Fail

Using a CH_3COOLi solution that contains 200 mM CH_3COOLi in 1% SDS, we were able to extract the gDNA from *Saccharomyces cerevisiae* efficiently. (Marko Lööke Article: Biotechniques.com Volume 50 No. 5 '11) "Biotechniques 50:325-328"

Protocol:

- 1) At log growth stage for your yeast colonies, use an inoculating loop to select a few colonies. (from a YPD growth plate)
- 2) Place the colonies selected into 100 μl 200 mM CH_3COOLi 1% SDS solution.
- 3) Incubate, after suspending cells, at 70°C for 15 minutes
- 4) Add 300 μl 95% ethanol
- 5) Vortex briefly to mix.
- 6) Centrifuge @ 15,000 g for 3 minutes
- 7) Dump supernatant.
- 8) Resuspend the precipitated DNA in 100 μl T.E. buffer.
- 9) Spin for 1 min @ 15,000 g (cellular debris will form a pellet at the bottom.)
- 10) Supernatant is contains gDNA from your yeast sample.
- 11) Load 1 μl into PCR tube for further testing/Amplification.

DNA Densities: * Jared: 11.1 ng/ μl /*
Chris: 13.1 ng/ μl
Neli: 12.3 ng/ μl

* We used for q-PCR
(Phusion)